



# GER-MA-WAY DISINFECTANT SPRAY

Hospital Grade, Ready to Use  
Clean, Disinfect and Protect Delicate  
Medical Devices

## วิธีการใช้งาน

น้ำยาฆ่าเชื้อพร้อมใช้สำหรับทำความสะอาดฆ่าเชื้อพื้นผิวอุปกรณ์ทางการแพทย์และเครื่องมือแพทย์ต่างๆพร้อมใช้ไม่ต้องผสมให้ประสิทธิภาพในการฆ่าเชื้อที่กว้างและรวดเร็ว มีผลทดสอบฆ่าเชื้อก่อโรคในโรงพยาบาลทั้งแบคทีเรียแกรมบวกและแกรมลบ ยีสต์และไวรัสรวมถึงโควิด-19 มี Contact time ไม่เกิน 1 นาที สามารถใช้ทำความสะอาดฆ่าเชื้อพื้นผิววัสดุอุปกรณ์เครื่องมือทางการแพทย์ได้หลากหลายโดยไม่ทำลายพื้นผิว Ger-Ma-Way ช่วยประหยัดเวลาและค่าใช้จ่ายในการควบคุมการติดเชื้อในสถานพยาบาลและสถานประกอบการที่ต้องมีการควบคุมการแพร่กระจายเชื้อ

Activity Claim	Test Protocol	Contact Time	Test Organisms active against
Bactericidal Claim 1-minute kill time	EN13727 - dirty conditions	1 minute	Ps. aeruginosa, E. coli, Staphylococcus aureus, Proteus vulgaris, Salmonella enterica, Enterococcus hirae, MRSA, VRE faecium, Klebsiella pneumoniae, A. baumannii, A. baumannii multi resistant drug
Fungicidal Claim 1-minute kill time	EN13624 - dirty conditions	1 minute	Trichophyton mentagrophytes
Yeasticidal Claim 1-minute kill time	EN13624 - dirty conditions	1 minute	Candida albicans
		1 minute	Candida auris
Mycobactericidal Claim (Tuberculocidal) 1-minute kill time	EN14348 - dirty conditions	1 minute	Mycobacterium terrae
Virucidal Claim 1-minute kill time	EN14476 - dirty conditions	1 minute	Poliovirus
		1 minute	Herpes virus
		1 minute	Adenovirus
		1 minute	Influenza A virus (H3N2), Influenza A (H1N1), Influenza B
		30 seconds	Human coronavirus for COVID-19

Product Code	Description and Size
100854601	GER-MA-WAY DISINFECTANT SPRAY 750 ML
100854746	GER-MA-WAY DISINFECTANT SPRAY 5 LITER

Efficacy testing from independent laboratory,  
data available upon request  
Material compatibility testing from independent laboratory,  
data available upon request